 **BodyMap**[®]
VACUUM SYSTEM BY ORTHODESIGN

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VACUUM POSITIONING SYSTEM BODYMAP FOR TRUNK, HEAD AND
LIMBS STABILIZATION

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
Models:

- Seat cushion **BodyMap A**
- U cut seat cushion **BodyMap A-**
- Sitting positioning cushion with thigh abductor wedge and lateral supports **BodyMap A+**
- Back cushion **BodyMap B**
- Backrest cushion with lateral supports **BodyMap B+**
- Back cushion with headrest and lateral support **BodyMap C**
- Positioning chair **BodyMap AB**
- Orthopedic seat for children with pelvic, trunk and head stabilization **BodyMap AB+**
- Positioning chair with headrest and lateral support **BodyMap AC**
- Headrest **BodyMap D**
- Neck cushion **BodyMap DW**
- Butterfly headrest **BodyMap DX**
- Straight headrest **BodyMap DY**
- Triangular headrest **BodyMap DZ**
- Upper limb support for wheelchairs, strollers and standing frames **BodyMap E**
- Lower limb support **BodyMap F**
- Derotation stabilizer **BodyMap G**
- Unilateral trunk stabilizer **BodyMap H**
- Head support **BodyMap I**
- Upper limb support **BodyMap J**
- Stabilizing mattress **BodyMap K**
- Backrest for car seat **BodyMap L**
- Cushion stabilizing the upper limbs **BodyMap M**
- Feet correction cushion **BodyMap N**
- Stabilizing roller **BodyMap O**
- Vacuum limb immobilizer **BodyMap P**
- Sensory bean bag floor cushion **BodyMap R**
- Sensory bean bag floor cushion **BodyMap R+**
- Bath ring **BodyMap S Swimmer**
- Frames **BodyMap**
- Vismemo thermoactive cover **BodyMap VisMemo**

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EN

BodyMap® System refers to first class medical products which are designed for adults as well as children with cerebral palsy. All these products provide perfect stabilization during seating or lying position. What is more, it is a perfect choice for the disabled due to the fact that it prevents bed sores.

BodyMap® System uses vacuum method. This method refers to products which are filled with polystyrene granules thanks to which products can model according to patient's body shape. This method is well-known within The Western European countries and it is widely applied in stabilization process. Furthermore, it is also applied during transportation of patients in ambulances. Polystyrene granules move freely inside the device and, more to the point, it is easy to move them in order to shape cushion at your discretion. The shape of cushion doesn't change easily because of the fact that the air is sucked out of the cushion.

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Concerning shaping the cushion, it is highly important to open valve so as to let the air in. This activity is done in order to equalize pressure inside the device with atmospheric pressure.

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
In order to open valve it is necessary to put thumb on frontal valve surface. Moreover, use your middle as well as index finger to press blue ring (1) and simultaneously slip off valve from pipe (2).

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
Be careful with the valve because it is small and it can be easily lost. Be cautious especially when you perform all activities outside.

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Cushion without valve should be left for few minutes only so as to equalize pressures.

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Next, shape cushion at your discretion. You can shape it by sitting or leaning on it - it depends on the type of cushion.

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
Once cushion is shaped it is necessary to let out the air from cushion.

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
If product was bought with an electric vacuum pipe it is necessary to connect these devices. Insert pipe's self-closing part into valve which is mounted on pipe.

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Large hand pump is added in case when buyer bought outsize product. Large hand pump connects with the product in the same way as small hand vacuum pump.

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
As far as standard products are concerned, manufacturer adds small hand pump to each product.

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EN

Move piston forward and backward so as to let out the air.

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
Air should be let out until it is properly hard. The shape of cushion can be easily shaped during this activity. Pump out longer so as to make cushion harder and more solid.

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
Once the shape as well as hardness of cushion is done, it is necessary to disconnect vacuum pump. Valve will automatically close air supply for cushion.

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Press metal valve metal ring backwards in order to disconnect pump. Next, take out vacuum pump's pipe.

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Cushion keeps its shape for a long time. If the shape of product changed, it is important to connect vacuum pump again so as to pump out the air. It is possible to buy thermo-active cover which makes the product more comfortable to use.

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BM_001 SMALL HAND VACUUM PUMP

Handy, small vacuum pump. Perfect for small and medium cushions. Each cushion includes this vacuum pump free of charge.

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BM_002 LARGE HAND PUMP

Very efficient manual vacuum pump. Perfect for mattresses and big cushions. Each mattress includes this vacuum pump free of charge.

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BM_003 ELECTRIC PUMP

Handy and very efficient electric vacuum pump. Equipped with electric charger. Can work with AA battery also (not included).

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BM_004 MANUAL DUAL-FUNCTION PUMP

Efficient, handy and made of high quality components. Strongly recommended for small cushions.

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BM_005 ELECTRIC DUAL-FUNCTION PUMP

With electric dual-function pump it is possible to adjust the hardness of the cushion to the individual needs of the patient. This pump allows precise pumping out or in the air in order to achieve the desired hardness. Noiseless, efficient and made of durable aluminium. Strongly recommended for large cushions and mattresses.

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EN

Cushions should be regularly cleaned with the aid of dry brush or vacuum cleaner equipped with attachment made of soft bristles. If stains are big, it is necessary to clean them with running water in a temperature up to 40 °C.

It is also possible to use detergent such as washing-up liquid or soap solution. Remember to rinse cleaned place off with running water and dry it with soft towel.

It is forbidden to use solvents: toluene or acetone. It is forbidden to use chlorine or its derivatives.

Dry far from direct heat.

Remember to be careful when you clean banderoles which protect places that connect cushion's elements.

It is necessary to be cautious not to cut or pierce cushion; do not use any sharp tools for cleaning.

You can repair damaged cushion using "repair kit" which is included free of charge to each product. Kit consist of special neoprene super glue and few patches which you can use to seal the damage surface of the cushion.

Cushion should be sent to manufacturer in case of mechanical damages.

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EN

ATTENTION: To avoid falling down the granules during molding the cushions with backrest (B, C, AB, AC, etc.) please follow the steps below:

- Place the cushion on flat surface
- Spread the granules inside the cushion evenly
- vacuum the cushion a little bit so it can still be floppy and moldable
- put the cushion on its place
- start the final molding seating on it or mold it manually

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EN

After pumping the air out, the cushion should keep its shape for at least 48 hours *. In case it loses its shape sooner, check the cushion yourself by carrying out **THE AIR LEAKAGE TEST**:

STEP 1. Check if the cushion has no visible mechanical damage, tearing or signs of wear and if the connective elements have not been damaged. If the cushion has a small tearing, you can try to fix it yourself using the “REPAIR KIT”, which is added to the cushion (follow the instructions included). If the damage is significant, contact the manufacturer to have the cushion repaired (warranty does not cover mechanical damage or tearing of the cushion).

STEP 2. Carry out the maintenance of the tube, which is used for pumping the air out:

- take the blue tube out of the cushion by pressing the blue clenching elements and remove the valve,
- cut off about 10 mm of the tube on both sides,
- insert the tube back into the cushion and put the valve on it.

After **STEP 2**, check again if the cushion is airtight by shaping the cushion and pumping the air out of it. If the cushion still does not keep the shape, proceed to **STEP 3**.

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EN

STEP 3. If the cushion is equipped with the service valve, check if the valve is airtight:

- open the valve cap with the use of a special key,
- check if the surface around the valve or the seal are not damaged,
- clean and wash the seal with soap or dishwashing liquid; put on and close the valve cap, paying particular attention that the granules do not get under the seal.

All of the abovementioned steps should be done carefully because the sudden or rapid movements may result in the granules pouring out of the cushion. Do not press the cushion when the valve cap is removed.

After STEP 3, check again if the cushion is airtight by shaping the cushion and pumping the air out of it. If the cushion still does not keep the shape, send the cushion back to the manufacturer to have it repaired.

**There are a lot of factors which affect the shape of the BodyMap® cushion, such as the kind of the cushion, its size, frequency of use, ambient temperature and the weight of the patient. The cushion should keep its shape from several days to several weeks. If the cushion is not used and does not keep its shape in the period of 48 hours, it might mean that it is damaged. To check that, carry out the air leakage test which is described above. If the test confirms that the cushion is damaged, inform the manufacturer about that fact in order to have the cushion repaired or replaced.*

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EN

BodyMap® valve is an innovative as well as easy solution used in **BodyMap®** system. It makes it possible to control the amount of granules in BodyMap cushions. It is highly advantageous because then it is possible to adjust cushion according to your needs.

Valve is situated at the back of cushion thanks to which it is not very visible and, moreover, it is not an obstacle while using the product. What is more, it is really easy to use it – you just need to unscrew cork and either pour in or take out specific amount of granules. Additional amount of granules is packed in airtight bags that are available in two sizes. Remember that the size of bag depends on the size of **BodyMap®** cushion.

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**EN**

Please follow these steps in order to pour in/take out granules:

1. Vacuum the cushion using the vacuum pump.

2. Unscrew carefully valve's cork.

CAUTION! When valve is unscrewed it is forbidden to move pillow violently or press it because then granules may spill out. There is a key added to the bag which facilitates opening and closing the service valve in order to add or decrease the amount of the granules

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3. Take out cork which secures funnel from bag with granules. Be careful so as not to spill out its contents.

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4. Insert tightly bag's funnel into hole from cushion at a depth of 1cm. Pull bag upwards, pour granules in as much as it is needed. If you want to take out some granules from cushion, it is necessary to insert funnel into hole from cushion and then turn cushion in a way that valve is placed down. Next, take out granules as much as it is needed.

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ATTENTION!

Before using the **BodyMap®** products, always make sure that the place of use (e.g. bed, mattress, seat) is free from sharp edges, metal parts, springs and other elements which may damage the product.

The manufacturer is not liable for the mechanical damage of the product!

ATTENTION!

In order to minimize the risk of bedsores, we recommend using a cotton cover along with each **BodyMap®** cushion if the permanent duration of use is longer than 2 hours. For people using **BodyMap®** cushions permanently for more than 8 hours it is recommended to use it along with thermo-elastic **VISmemo™** cover. Caution should be taken in subjects to people with known bedsores - prone. In such a case it is strictly recommended to you use **BodyMap®** cushion with aeration-massaging pad for people who use **BodyMap** cushions non-stop longer than 8 hours!

Under any circumstances BodyMap® product cannot be fold, when the air is pump out and the vacuum valve is assembled.

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ATTENTION!

Over time, the polystyrene granules, which fill the **BodyMap®** cushions, may give signs of flattening. As a result, the volume of the cushions will decrease, which is their natural feature. If this situation arises, we recommend replacing or refilling the granules in the cushion. It can be easily done with a use of a service valve which most **BodyMap®** cushions are equipped with. In particular, all seats, mattresses and beans.

In the cushions not equipped with a service valve it is not possible to refill the granules. We recommend you to consider the option of adding such a valve while ordering. The decreasing volume of the cushions is their natural feature and it shall not give rise to a complaint.

Detailed information about the **BodyMap®** cushions equipped with a service valve can be found in the catalogue and on the website www.bodymapsystem.com

ATTENTION!

First, place the **BodyMap®** cushion on the chair, stroller etc.

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ATTENTION!

It is not allowed to press and load the **BodyMap®** device in case the air is inside and the vacuum valve is assembled. In order to shape the cushion, there is necessary to spread the granulate evenly and to pump the air out slightly. By maintaining the cushion soft, you can shape the device according to your needs.

ATTENTION!

The Velcro surface of **BodyMap®** cushions is used for positioning of cushions on the product or fix the supporting elements on it.

ATTENTION!

There are a lot of factors which affect the shape of the **BodyMap®** cushion, such as the kind of the cushion, its size, frequency of use, ambient temperature and the weight of the patient. The cushion should keep its shape from several days to several weeks. If the cushion is not used and does not keep its shape in the period of 48 hours, it might mean that it is damaged and the manufacturer should be informed about that fact in order to repair or replace the faulty cushion.

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